



PRODUCT PORTFOLIO

RIGID GAS PERMEABLE CONTACT LENSES

SOFT CONTACT LENSES

LENSCARE SOLUTIONS

ACCESSORIES

MENICON SAS

MENICON SERVICES

■ MENICON WEBSITES

www.menicon.com • www.miru.menicon.com • www.solocareaqu.eu

■ CUSTOMER SERVICE

From Monday to Thursday 9:00 AM - 6:00 PM
Friday 9:00 AM - 5:00 PM

Tel.: +33 1 41 40 19 19

Fax: +33 1 41 40 19 20
serviceclients@menicon.fr

■ TECHNICAL ASSISTANCE

From Monday to Thursday 9:00 AM - 6:00 PM
Friday 9:00 AM - 5:00 PM

Tel.: +33 1 41 40 19 41

Fax: +33 1 41 40 19 20
atechnique@menicon.fr

■ TRAINING AND MARKETING TOOLS

Contact your Menicon regional sales manager

For a new product, please inquire with your Menicon Regional Sales Manager or Menicon Customer Service about the conditions. Regulatory confirmation may be required before taking any order by Menicon.



Infinite Vision


From the time we open our eyes in the morning
until we fall asleep at night,
We use our eyes to see the world.
We create, we learn, we find new friends.
We take vision for granted,
but it enriches every day of our lives.
Here at Menicon, we contribute to daily living
by providing superior visual correction:
Working to bring each and every person the joy of sight.



SOFT CONTACT LENSES

Replacement		1 Day	\emptyset_T (mm)	r_0 (mm)	Power (D)	Cylinder (D) Axis (°)
<p>Miru 1day Menicon Flat Pack</p> <p>Box of 30 lenses</p> <p>Material: Hioxifilcon A (HEMA - GMA) 57 % water content Dk/e : 25×10^{-9}</p>		14.20	8.60	+4.00 to -6.00 D in 0.25 D steps -6.50 to -10.00 D	-	

Replacement		1 Month	\emptyset_T (mm)	r_0 (mm)	Power (D)	Cylinder (D) Axis (°)
<p>Miru 1month Menicon</p> <p>Box of 6 lenses</p> <p>Material: Silicone hydrogel Asmofilcon A 40% water content - Dk 129×10^{-11} Dk/e : 161×10^{-9}</p>		14.00	8.30 8.60	+6.00 to -6.00 D in 0.25 D steps -6.50 to -13.00 D in 0.50 D steps	-	
<p>Miru 1month Menicon For Astigmatism</p> <p>Box of 6 lenses</p> <p>Material: Silicone hydrogel Asmofilcon A 40% water content - Dk 129×10^{-11} Dk/e : 161×10^{-9}</p>		14.00	8.60	00.00 to -6.00 D in 0.25 D steps -6.50 to -10.00 D in 0.50 D steps	Cylinders -0.75 D -1.25 D -1.75 D Axis 10°, 20°, 90° 160°, 170°, 180° (80°, 100°: available in 2Q 2016)	

Replacement		1 Month	\emptyset_T (mm)	r_0 (mm)	Power (D)	Cylinder (D) Axis (°)
<p>Menisoft</p> <p>Box of 3 lenses</p> <p>Material: Hydrogel Mipafilcon A 72% water content - Dk 35×10^{-11} Dk/e : 43×10^{-9}</p>		14.20	8.60	-0.25 to -6.00 in 0.25 D steps -6.50 to -10.00 in 0.50 D steps	-	


Menicon recommends to use SOLOCARE AQUA® or MeniCare Soft™, multi-purpose solution


SOFT CONTACT LENSES

Replacement **2**
Weeks




		\varnothing_T (mm)	r_0 (mm)	Power (D)	Cylinder (D) Axis (°)
Menicon PremiO Box of 6 lenses Material: Silicone hydrogel Asmofilcon A 40% water content - Dk 129 x10 ⁻¹¹ Dk/e : 161x10 ⁻⁹		14.00	8.30 8.60	+6.00 to -6.25 D in 0.25 D steps	-
				-6.50 to -13.00 D in 0.50 D steps	
Menicon PremiO Toric Box of 6 lenses Material: Silicone hydrogel Asmofilcon A 40% water content - Dk 129 x10 ⁻¹¹ Dk/e : 161x10 ⁻⁹		14.00	8.60	Plan to -6.00 D in 0.25 D steps	Cylinders -0.75 D -1.25 D -1.75 D Axis 10°, 20°, 90° (except -1,75 D), 160°, 170°, 180°
				-6.50 to -10.00 D in 0.50 D steps	

Replacement **3**
Months

		\varnothing_T (mm),	r_0 (mm)	Power (D)	Cylinder (D) Axis (°)
Range Individual 1 lens or a 1 box of 2 lenses Individual Individual Toric Individual Progressive 2* (Central Distance) Individual Progressive 2 Toric Individual Progressive 2 Near (Central near) Individual Progressive 2 Near Toric		13.00 to 15.50 per 0.10 mm	7.40 to 9.50 per 0.10 mm	-25.00 to +25.00 D in 0.25 D steps Addition +1.00 to +3.00 in 0.50 D steps	Cylinders -0.75 D to -5.00 D -0.75 D to -3.00 D (Prog. 2 Toric and Prog. 2 Near Toric) in 0.25 D steps Axis from 0° to 180° per 1°
				Material: Silicone hydrogel Filcon V3 72% water content - Dk 60 x10 ⁻¹¹ *Individual Progressive is also available.	

Keratoconus contact lenses and special cases		
ROSE K2 Soft™ 1 lens or 1 box of 2 lenses with same parameters		(Ref page 9)

Replacement **1**
Year

	\varnothing_T (mm)	r_0 (mm) per 0.30 mm		Power (D) per 0.25 (D)	
		Standard	Non-standard	Standard	Non-standard
Menicon soft 72 Material: Mipafilcon A 72 % water content 	13.50	8.40 / 8.70 / 9.00	8.10	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00
	14.00	8.40 / 8.70 / 9.00	9.30		
Menicon soft 72 Toric Material: Mipafilcon A 72 % water content 	14.00	8.40 / 8.70 / 9.00		-1.00 to -8.00 Axis 10° to 180°* per 10°	
	Cyl : -0.75 to -2.75 in 0.50 D steps *Except: -1.00 and -1.25 D : 10° 20° 160° 170° 180°				
Menicon soft MA Material: Govafilcon A 38 % water content 	13.00	8.10 to 9.00	7.80 / 9.30 / 9.60	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00
	13.50	8.10 to 9.00	9.30 / 9.60 / 9.90		
	14.00	8.40 / 8.70 / 9.00	9.30 à 9.90		

Menicon recommends to use SOLOCARE AQUA® or MeniCare Soft™, multi-purpose solution

LENSCARE SOLUTIONS

FOR ALL SOFT CONTACT LENSES

■ SOLOCARE AQUA® multi-purpose solution with MicroBlock® lens case



Tri Pack SOLOCARE AQUA®

3 x 360ml + 3 lens case
Order unit: 6



Duo Pack SOLOCARE AQUA®

2 x 360ml + 2 lens case
Order unit: 12



SOLOCARE AQUA®

360ml + lens case
Order unit: 10



SOLOCARE AQUA®

90ml + lens case
Order unit: 12



Travel kit SOLOCARE AQUA®

90ml + lens case + mirror
Order unit: 12

■ MENICARE SOFT multi-purpose solution



Pack MeniCare Soft

2 x 360ml + 2 lens case
Order unit: 10



MeniCare Soft

360ml + lens case
50ml + lens case / 70ml + lens case*
Order unit: 10

* the volume of small bottle will be changed from 50ml to 70ml

ACCESSORIES

FOR ALL SOFT CONTACT LENSES

HEMECASE Miru 1day Menicon Flat Pack

Storage unit
30 "Flat Pack" (R) / 30 "Flat Pack" (L) can be stored.



Meniholder

Order unit 10 holders



RIGID GAS PERMEABLE CONTACT LENSES

SPHERIC CONTACT LENSES

	\varnothing_T (mm)	r_0 (mm) in 0.05 mm steps		Power (D) in 0.25 (D) steps		
		Standard	Non-standard	Standard	Non-standard	
Z MATERIAL	Menicon Z Aspheric Material : Tisilfocon A Dk: 163 (ISO)	9.00 9.20 9.40 9.60 9.80	7.20 7.30 to 7.50 7.60 to 7.80 7.90 to 8.10 8.20 to 8.40 in 0.10 mm steps	\varnothing_T from 9.00 to 9.80 for 7.20 to 8.60 in 0.10 mm steps	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00
	Menicon Z-α Sphero-aspheric	8.00 to 9.00 in 0,20 mm steps 8.80 9.20 9.40 9.60 9.80 10.00 10.50 11.00	7.20 to 8.40 in 0.05 mm steps	6.50 to 9.00 in 0.05 mm steps 6.70 to 9.00 7.30 to 9.00 7.50 to 9.00	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00
	Menicon EX-Z Sphero-aspheric	8.80 9.20 9.60 10.00 10.50 11.00*	7.20 to 8.40 in 0.05 mm steps	6.50 to 8.60 6.50 to 8.80 6.50 to 9.00 6.90 to 9.00 7.20 to 9.00	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00
	Menicon Z Comfort with tear-reservoir zone	10.20 10.60 11.00	7.20 to 8.40 in 0.05 mm steps	6.50 to 9.00 in 0.05 mm steps	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00
	Menicon Z OMNI Sphero-aspheric (Z or EX material)	9.10 9.60 10.10	7.20 to 8.40 mm	6.50 to 9.00 in 0.05 mm steps	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00
EX MATERIAL	Menicon EX Sphero-aspheric Material: Tolofocon A DK: 64 (ISO)	8.80 9.20 9.60 10.00 10.50 11.00*	7.20 to 8.40 in 0.05 mm steps	6.50 to 8.60 6.50 to 8.80 6.50 to 9.00 6.90 to 9.00 7.20 to 9.00	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00

* Lenticular geometry EX 11.00 mm \varnothing_0 : 7.00 mm

TORIC CONTACT LENSES

	\varnothing_T (mm)	r_0 (mm) in 0.05 mm steps		Power (D) in 0.25 (D) steps		
		Standard	Non-standard	Standard	Non-standard	
Z OR EX MATERIAL	Menicon Z BTC (Back toric compensation)	9.10 9.60 10.10	7.20 to 8.40 mm	7.00 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00
	Menicon Z Comfort BTC (Back toric with large diameter and tear-reservoir zone)	10.20 10.60 11.00	7.20 to 8.40 mm	6.50 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00
	Menicon Z BT (Bitoric)	9.10 9.60 10.10	7.20 to 8.40 mm	7.00 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00
	Menicon Z Comfort BT (Bitoric with large diameter and tear-reservoir zone)	10.20 10.60 11.00	7.20 to 8.40 mm	7.00 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00
	Menicon Z FT (Front Toric)	9.10 9.60 10.10	7.20 to 8.40 mm	6.50 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00

* Other parameters might be available upon request, please contact our Technical Assistance.

Menicon recommends to use MENICARE PURE or MeniCare Plus, multi-purpose solution for daily care + PROGENT once a week.

RIGID GAS PERMEABLE CONTACT LENSES

MULTIFOCAL CONTACT LENSES

ZOR EX MATERIAL		\varnothing_T (mm)	r_0 (mm) per 0.05 mm		Power (D) per 0.25 (D)		
			Standard	Non-standard	Standard	Non-standard	
SPHERIC	Simultaneous vision						
	Menicon Z Progressive (Center Distance)	9.10	7.20 to 8.40 mm <i>Addition : +1.00 to +3.00 in 0.50 D steps</i>	6.50 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00	
		9.60					
		10.10					
	Z Comfort Progressive (Center Distance, large diameter, tear-reservoir zone)	10.20	7.20 to 8.40 mm	6.50 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00	
		10.60					
		11.00					
	Z Comfort Progressive Near	10.20	<i>Addition : +1.00 to +3.00 in 0.50 D steps</i>				
		10.60					
		11.00					
	Z Progressive Near (Center Near)	9.10	7.20 to 8.40 mm	6.50 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00	
		9.60					
		10.10					
	Z Progressive + (accelerated near vision progression)	9.10	<i>Addition : +1.00 to +3.00 in 0.50 D steps</i>				
		9.60					
10.10							
BACK TORIC	Z Progressive BTC	9.10	7.20 to 8.40 mm	7.00 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00	
		9.60	<i>Addition : +1.00 to +3.00 in 0.50 D steps</i>				
		10.10	<i>Back cylinder (mm) : 0.20 to 1.50 in 0.05 steps</i>				
	Z Comfort Progressive BTC	10.20	7.20 to 8.40 mm	6.50 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00	
		10.60	<i>Addition : +1.00 to +3.00 in 0.50 D steps</i>				
		11.00	<i>Back cylinder (mm) : 0.20 to 1.00 in 0.05 steps</i>				
	Z Progressive + BTC	9.10	7.20 to 8.40 mm	7.00 to 9.00 mm	+5.00 to -10.00		
		9.60	<i>Addition : +1.00 to +3.00 in 0.50 D steps</i>				
		10.10	<i>Back cylinder (mm) : 0.20 to 1.50 in 0.05 steps</i>				
BITORIC	Z Progressive BT	9.10	7.20 to 8.40 mm	7.00 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00	
		9.60					
		10.10					
	Z Comfort Progressive BT	10.20	7.20 to 8.40 mm	7.00 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00	
		10.60					
		11.00					
Z Progressive + BT	9.10	<i>Addition : +1.00 to +3.00 in 0.50 D steps</i>					
	9.60	<i>Back cylinder (mm) : 0.20 to 1.50 in 0.05 steps (except for Z Comfort Prog. BT : 0.20 to 1.00)</i>					
	10.10	<i>Front cylinder (D) : -0.50 to -2.00 in 0.25 steps, Axis: 0° to 180° in 1° steps</i>					
	Alternating vision						
	Menifocal Z (Z material only)	9.60	7.20 to 8.40 mm per 0.05 mm	7.00 to 9.00 mm	+5.00 to -10.00	-10.25 to -13.00	
		9.80					<i>Addition : +1.00 to +2.50 in 0.50 D steps</i>
Menicon Z Executive (segmented)	9.10	7.20 to 8.40 mm per 0.05 mm	6.50 to 9.00 mm	+5.00 to -10.00	-10.25 to -25.00		
	9.60					<i>Addition : +1.00 to +3.00 in 0.50 D steps</i>	
	10.10					<i>Prisme ballast : 1.00 to 2.50 in 0.50 Δ steps (standard : 2.00) Position of segment : +0.5mm / -0.5mm</i>	

Upon request, lenses are available with Dk 58 material (except for Menifocal Z).

* Other parameters might be available upon request, please contact our Technical Assistance.

ORTHOKERATOLOGY (Corneal reshaping overnight lens)

Z MATERIAL	\varnothing_T (mm)	r_0 (mm) 7.00 to 10.00 mm per 0,05 mm	
Menicon Z Night	10.20	<i>Tangential angle: 46° to 63° (in 1° steps)</i> <i>Sagittal depth: 0.30 to 1.00 mm (in 0.01mm steps)</i>	
	10.60		
Menicon Z Night Toric	11.00	<i>Tangential angle: 46° to 63° (in 1° steps)</i> <i>In the flat meridian (FL) and in the steep meridian (ST)</i>	
		<i>Sagittal depth: 0.30 to 1.00mm (in 0.01mm steps)</i> <i>In the flat meridian (FL) and in the steep meridian (ST)</i>	

Menicon recommends to use MENICARE PURE or MeniCare Plus, multi-purpose solution for daily care + PROGENT once a week.

RIGID GAS PERMEABLE LENSES - OTHER DESIGNS

	\varnothing_T (mm)	r_0 (mm) in 0.05 steps	F'v per 0.25 (D)	
Z O R E X M A T E R I A L	Menicon Plateau	8.00 to 11.00 in 0.10 steps (standard : 10.00)	7.03 to 10.23 (48.00 to 33.00 D) r_1 1.00 to 7.00 D per 1.00 D	+20.00 to -20.00
	Menicon S9	8.50 8.80 9.00 9.20 9.50	6.50 to 9.00	+25.00 to -25.00
	Menicon KRC 1 to 12	7.90 8.00 8.10	5.15 to 6.90	-0.25 to -25.00
	Menicon K3	8.50	5.00 to 6.90	-0.25 to -25.00
	Menicon TCC6	9.20 9.30 9.50 10.00 10.50	6.50 to 9.00	+25.00 to -25.00
	Menicon B4P	10.00 10.50 11.00 11.50 12.00	6.50 to 9.00 7.20 to 9.00 7.25 to 9.00	+25.00 to -25.00
	Menicon B4PM	11.00 11.50 12.00	7.20 to 9.00 7.25 to 9.00	+25.00 to -25.00
	Menicon MX	8.50 8.80 9.00	6.50 to 9.00	+25.00 to -25.00
Other designs, please contact our Technical Assistance.				

Menicon recommends to use MENICARE PURE or MeniCare Plus, multi-purpose solution for daily care + PROGENT once a week.

SOFTWARE FOR EYE CARE PROFESSIONALS

easyfit™

Fit Gas Permeable & Soft Contact Lenses
quickly & easily


To get more information: Contact your Menicon regional sales manager

Three functions:

- 1 Recommend suitable lenses based on Topography data
 - RGP spherical, toric, multifocal
 - ROSE K2 Range
 - Menicon ZNight, ZNight Toric
 - Miru 1day Menicon Flat Pack
 - Menicon PremiO, PremiO Toric
 - Miru 1month Menicon
 - Miru 1month Menicon For Astigmatism
 - Individual Range
- 2 Order directly from easyfit
- 3 Get technical advice more easily
 - Send the patient data file just with one click

Compatible with more than 25 topographers

BON Eyetop	Phoenix
Cassini	Rodenstock DNEye Scanner
EyeTop 2005 Haag-Streit	Shin-Nippon CT-1000
Keratron	Tomey RT-7000
Essilor Anaeyes	Tomey TMS-2 / TMS-3 / TMS-4
Medmont (version 3)	Topcon 7000
Medmont (version 4)	Topcon 8100 (CA-100 software)
Medmont (version 5)	Topcon 8100 (Color Mapping software)
Medmont (version 6)	Topcon CA-100 / CA-200 / CA-800
Nidek OPD-Scan III	Visionix VX-120
Oculus Easygraph	Zeiss i-Profiler
Oculus Keratograph	Zeiss i-Profiler Plus
Oculus Pentacam	



Realistic and dynamic fluorescein images

With Menicon design lenses

Access to technical documentations

Information: www.e-menicon.fr

KERATOCONUS IRREGULAR CORNEA

ROSE K LENSES



R O S E K

	\varnothing_T (mm)	r_0 (mm)	F'v per 0.25 (D)	Edge Lift (mm)	
Z O R E X MATERIAL	ROSE K2 (Aberration control) ROSE K	7.90 to 10.40 in 0.10 steps (standard : 8.70)	4.30 to 8.55 in 0.05 steps ¹ from +3 decreased to -1.3 increased in 0.10 increments	+ 25.00 to -25.00 standard 0.0 flat +1.0 steep -0.5 others ¹	
	ROSE K2 POST GRAFT	9.40 to 12.00 in 0.10 steps (standard : 10.40)	5.70 to 9.30 in 0.10 steps	+25.00 to -25.00 standard 0.0 flat +0.6 steep -0.5 double flat +1.2 double steep -1.0	
	ROSE K2 IC	9.40 to 12.00 in 0.10 steps (standard : 11.20) 11.20, $r_0 \geq 6.00$ 12.00, $r_0 \geq 6.70$	5.70 to 9.30 in 0.10 steps	+25.00 to -25.00 standard 0.0 flat +0.6 steep -0.5 double flat +1.2 double steep -1.0	
	ROSE K2 NC	7.60 to 9.00 in 0.10 steps (standard : 8.30)	4.30 to 7.60 in 0.05 steps	+10.00 to -35.00 in 0.25 steps ² Available variations -1.5 to +3.00 in 0.1 steps	standard 0.0 steep -0.5 flat +1.0 others ²
	ROSE K2 XL (Z material only)	13.60 to 15.60 in 0.10 steps (standard : 14.60)	5.80 to 8.40 in 0.10 steps	+15.00 to -30.00 in 0,25 steps ³ from +3.00 decreased to -3.00 increased in 0.50 increments	standard 0.0 flat +1.00 steep -1.00 double flat +2.0 double steep -2,00 others ³

ROSE K2, ROSE K2 Post Graft, ROSE K2 IC
ROSE K2 XL) Option toric or asymmetric:

- Toric peripheral curves: 0.40 to 2.00mm.
- Asymmetric corneal technology (ACT): 0.40 to 1.50mm (in 0.10mm steps).
(0.80 to 2.00mm for ROSE K2 XL).
- Toric back, front and bi-toric surfaces, axis : 0° to 180° (in 1° steps).

For material availability, other possibilities of parameters, asymmetric design, please contact our Technical Assistance:
Phone +33 1 41 40 19 41 - atehnique@menicon.fr - ROSE K2 XL: Z Material only.

Menicon recommends to use MENICARE PURE or MeniCare Plus, multi-purpose solution for daily care + PROGENT once a week.

SOFT LENS FOR IRREGULAR CORNEAS

ROSE K2 Soft is a daily wear soft lens for irregular corneas, 3 month replacement lens manufactured from Silicon Hydrogel material, or 6 month or annual replacement lens manufactured from Hydrogel material.



R O S E K

	\varnothing_T (mm)	e_C (mm)	r_0 (mm)	F'v in 0.25 (D) steps	Edge Lift (mm)
ROSE K2 Soft™ Materials : - Silicone hydrogel Lagado (LSiH), Dk 49 49% water content - Menicon soft 72 Hydrogel Dk 34 72% water content	14.30 to 15.30 in 0.10 steps (standard : 14.80)	0.25 to 0.60 in 0.10 steps (standard : 0.35)	7.40 to 9.00 in 0.20 steps	-30.00 to +30.00 Cylinders -0.25 D to -10.00 D in 0.25 D steps Axis from 0° to 180° per 1°	standard 0.0 flat +1.0 steep -1.0 double flat +2.0 double steep -2.0

Menicon recommends to use Menicare Soft™ or SOLOCARE AQUA®, multi-purpose solution

LENSCARE SOLUTIONS

FOR ALL RIGID GAS PERMEABLE CONTACT LENSES

■ MENICARE PURE (multi-purpose solution)

MeniCare Pure

Multi-purpose solution
(daily use)

250ml + lens case

70 ml + lens case

Order unit: 10



■ PROGENT (protein remover)

New packaging

Menicon Progent (5 treatments)

Intensive cleaner for protein and for disinfection
(weekly use)

Box of 5 doses A + 5 dose B + 1 Progent vial

Order unit: 10



■ NEW INSEPARABLE PACKS (MENICARE PURE & PROGENT)

Starter Kit

MeniCare Pure & Progent

MeniCare Pure 70ml 1 bottle

+ 1 lens case, 1 dose of Progent A&B

+ 1 Progent vial

Order unit: 10



MENICARE PURE

MULTI-PURPOSE SOLUTION

for all rigid gas permeable lenses


Indications

- Cleaning
- Disinfecting: effective against bacteria and fungi in 30 minutes and active against acantamoeba in 4 hours
- Rinsing
- Storage

Composition

ε -Polylysine 0.01%
Edetate disodium 0.05%
HCO 60 0.1%

NATURAL FORMULATION



LENSCARE SOLUTIONS

FOR ALL RIGID GAS PERMEABLE CONTACT LENSES

■ MENICARE PLUS (multi-purpose solution)

MeniCare Plus

Multi-purpose solution (daily use)

250ml + lens case

Order units: 10



■ SPRAY & CLEAN (cleaning solution)

Spray & Clean

Cleaning solution for all rigid gas permeable lenses
Preservative free solution 15ml
Anti-lipid action

Order unit: 10



■ MENILAB (for trial lenses)

Menilab

For professional use
Disinfection of RGP trial lenses
250ml

Order unit: 1



ACCESSORIES

FOR ALL RIGID GAS PERMEABLE CONTACT LENSES

RGP lens case (Menicase)

individual lens case

Order unit: 10



Suction cup (Menicup)

Order unit: 10



Travel Kit (Menibox)

Case with mirror,
2 lens case, 1 suction cup

Order unit: 1



FOR SEMI-SCLERAL CONTACT LENSES (ROSE K2 XL)

Suction cup DMV® Maxi-Cup™

(To put a lens)

Order unit: 3



Suction cup DMV® Ultima™

(To remove a lens)

Order unit: 10



Large diameter RGP lens case (PROGENT treatment)

(PROGENT treatment)

Order unit: 1



Flat lens case

(For large diameter lenses, storage case)

Order unit: 1





www.menicon.com

MENICON SAS - 104 RUE MARTRE - 92583 CLICHY CEDEX - FRANCE
S.A.S AU CAPITAL DE 4 023 900 EUROS - RCS NANTERRE B 722 038 668 - CODE INTRACOMMUNAUTAIRE FR92722038668