



Simply find the solution

From the wide range of contact lenses

PRODUCT PORTFOLIO

RIGID GAS PERMEABLE CONTACT LENSES,
SOFT CONTACT LENSES, LENS CARE SOLUTIONS, ACCESSORIES

Summary:

Soft Contact Lenses	3 to 7	Lens Care Solutions & Accessories	12 to 14
Rigid Gas Permeable Contact Lenses	8 to 9	easyfit , topographer	15
Orthokeratology, Scleral Contact Lenses	10	Services	16
Keratoconus & Irregular Cornea	11		


Simply find the solution

From the wide range of contact lenses

All types of lenses	Soft Lenses	RGP Lenses
All modalities of replacement	1 day, 2 weeks, 1 month, 3 months	1 year
All types of wearing	Occasional, Daily wear, Extended wear (PremiO)	Daily to extended wear (Menicon Z)
All indications	Myopia, Hyperopia, Astigmatism, Presbyopia, Irregular cornea, Keratoconus	Myopia, Hyperopia, Astigmatism, Presbyopia, Irregular cornea, Keratoconus
Associated lens care solutions	MeniCare Soft, SOLOCARE AQUA®	MENICARE PURE or MENICARE PLUS PROGENT

All lenses, lens care solutions and accessories manufactured and / or distributed under the brand name Menicon or ROSE K are medical devices class IIa, IIb, I, respectively, marked CE (CE 0483, CE 0344, CE 1639). Menicon SAS is certified with NF EN ISO 13485.

TO ORDER ONLINE MENICON PRODUCTS



E-menicon website only for Eye Care Professionals:

- Order your Menicon products anytime
- Stay tuned about latest news

www.e-menicon.fr

Miru

1day Menicon Flat Pack



Replacement

1
Day

Material: Hydrogel Hioxifilcon A - Water content: 57% - Dk/t : 19 x 10 ⁻⁹ (ISO)		Ø _T (mm)	BC (mm)	Power (D)	Cylinder (D) Axis (°)
Miru 1day Menicon Flat Pack Daily wear 30-lens pack (daily disposable lens)		14.20	8.60	+4.00 to +0.50 per 0.25 D -0.50 to -6.00 per 0.25 D -6.50 to -10.00 per 0.50 D	
Miru 1day Menicon Flat Pack Daily wear 6-lens pack (daily disposable lens)					
Miru 1day Menicon Flat Pack toric Daily wear 30-lens pack (daily disposable lens)		14.50	8.60	Plano to -6.00 per 0.25 D*	Cylinders: -0.75 D -1.25 D -1.75 D** Axis: 15°, 90°, 165°, 180°
Miru 1day Menicon Flat Pack toric Daily wear 6-lens pack (daily disposable lens)					

New

New

ABOUT

Miru

1day Menicon Flat Pack

Daily disposable contact lenses
with Smart Touch™ and
Smart Fit technology

The slimmest contact lens packaging...


...in the world



- ✓ **Easy to handle**
- ✓ **Less raw material**
- ✓ **Less waste**
- ✓ **Recyclable blister**

New

Miru 1day Menicon Flat Pack toric



- **Auto-orientation** toric lens: Smart Fit technology
- Minimal average thickness profile
- Anatomically matches eyelids
- High-level comfort and axis stability

Smart Fit technology

The lens naturally orients in
the right way,
no matter in
which way it
is inserted.



Miru

1day UpSide

Replacement

1
Day

Material: Silicone Hydrogel Midafilcon A - Water content: 56% - UV Filter - Dk/t 91 x 10 ⁻⁹ (ISO) Ultra low modulus 0,36 Mpa - Low friction coefficient (0,023)		Ø _T (mm)	BC (mm)	Power (D)
Miru 1day UpSide Daily wear 30-lens pack (daily disposable lens)		14,20	8,40	-0,25 to -6,00 D per 0,25 D -6,50 to -10,00 D per 0,50 D
Miru 1day UpSide Daily wear 5-lens pack (daily disposable lens)				
Miru 1day UpSide Daily wear 3x30-lens pack (daily disposable lens)				
Miru 1day UpSide Daily wear 5x30-lens pack (daily disposable lens)				
Miru 1day UpSide Daily wear 10 x 30 lens R & L BOX "COLLECTION"				

SMART TOUCH™

Daily disposable contact lenses always in the right way
and ready to wear

Lens case with Smart Touch™
technology

type.1 Miru 1day UpSide




type.2 Miru 1day Menicon
Flat Pack



Menicon Smart Touch™

SMART TOUCH



<https://www.menicon.com/smarttouch/global/>



- Hygienic
- Hassle free
- Quick handling

Smart Touch™:
Healthy clean lenses
Hassle free handling!

Health - Vision
Comfort - Convenience

Replacement

2
Weeks

Material: Silicone Hydrogel Asmofilcon A - Water content: 40% Dk/t 161 x 10 ⁻⁹ (ISO)		Ø _T (mm)	BC (mm)	Power (D)	Cylinder (D) Axis (°)
Menicon PremiO Daily or extended wear (Up to 7 days / 6 nights) 6-lens pack		14.00	8.30 8.60	+6.00 to +0.25 per 0.25 D -0.25 to -6.25 per 0.25 D -6.50 to -13.00 per 0.50 D	-
Menicon PremiO Toric Daily or extended wear (Up to 7 days / 6 nights) 6-lens pack		14.00	8.60	Plano to -6.00 per 0.25 D -6.50 to -10.00 D per 0.50 D	Cylinders: -0.75 D -1.25 D -1.75 D Axis: 10°, 20°, 90° (except -1.75 D), 160°, 170°, 180°

Menicon recommends to use MeniCare Soft or SOLOCARE AQUA® multipurpose solution.

Miru 1month :

Complete range
of silicone hydrogel
monthly contact lenses



Replacement

1
Day

Material: Silicone Hydrogel Asmofilcon A - Water content: 40% Dk/t 161 x 10 ⁻⁹ (ISO)		Ø _T (mm)	BC (mm)	Power (D)	Cylinder (D) Axis (°)
Miru 1month Menicon Daily wear 6-lens pack		14.00	8.30 8.60	+6.00 to +0.25 per 0.25 D -0.25 to -6.00 per 0.25 D -6.50 to -13.00 per 0.50 D	-
Miru 1month Menicon toric Daily wear 6-lens pack		14.00	8.60	+4.00 to +0.25 per 0.25 Plano to -6.00 D per 0.25 D -6.50 to -10.00 D per 0.50 D	Cylinders: -0.75 D -1.25 D -1.75 D -2.25 D* Axis : 0° to 180° per 10° *-2.25 D only available in 10°, 20°, 90°, 160°, 170° and 180°
Miru 1month Menicon multifocal Daily wear 6-lens pack		14.20	8.60	+6.00 to +0.25 D per 0.25 D Plano to -6.00 D per 0.25 D -6.50 to -13.00 D per 0.50 D	LOW Design: Add. ≥ +1.00 D HIGH Design: Add. > +2.00 D

Menicon recommends to use MeniCare Soft or SOLOCARE AQUA® multi-purpose solution

SOFT CONTACT LENSES



Replacement

3
Months

Made to order lenses Material: Silicone hydrogel Filcon V3 74% water content - UV Filter - Dk 60 x 10 ⁻¹¹ (ISO)	Ø _T (mm)	BC (mm)	Power (D)	Cylinder (D) Axis (°)
Individual Daily wear 1 lens or a 2-lens pack	13.00 to 15.50 per 0.10 mm	7.40 to 9.50 per 0.10 mm	-25.00 to +25.00 D per 0.25 D	-
Individual Toric Daily wear 1 lens or a 2-lens pack	13.00 to 15.50 per 0.10 mm	7.40 to 9.50 per 0.10 mm	-25.00 to +25.00 D per 0.25 D	-0.75 D to -5.00 D per 0.25 D Axis: 0° to 180° per 1°
Individual Progressive 2* (Central Distance) Daily wear 1 lens or a 2-lens pack	13.00 to 15.50 per 0.10 mm	7.40 to 9.50 per 0.10 mm	-25.00 to +25.00 D per 0.25 D Addition +1.00 to +3.00 per 0.50 D	-
Individual Progressive 2 Toric (Central Distance) Daily wear 1 lens or a 2-lens pack	13.00 to 15.50 per 0.10 mm	7.40 to 9.50 per 0.10 mm	-25.00 to +25.00 D per 0.25 D Addition +1.00 to +3.00 per 0.50 D	-0.75 D to -3.00 D per 0.25 D Axis: 0° to 180° per 1°
Individual Progressive 2 Near (Central near) Daily wear 1 lens or a 2-lens pack	13.00 to 15.50 per 0.10 mm	7.40 to 9.50 per 0.10 mm	-25.00 to +25.00 D per 0.25 D Addition +1.00 to +3.00 per 0.50 D	-
Individual Progressive 2 Near Toric (Central near) Daily wear 1 lens or a 2-lens pack	13.00 to 15.50 per 0.10 mm	7.40 to 9.50 per 0.10 mm	-25.00 to +25.00 D per 0.25 D Addition +1.00 to +3.00 per 0.50 D	-0.75 D to -3.00 D per 0.25 D Axis: 0° to 180° per 1°

Menicon recommends to use MeniCare Soft or SOLOCARE AQUA® multipurpose solution.

*Individual Progressive is also available.

ROSE K2 SOFT™

Replacement

3
Months

6
Months

Rose K2 Soft is a **daily wear** soft lens for irregular corneas, 3 months replacement lens manufactured from Silicone Hydrogel material, or 6 months or annual replacement lens manufactured from Hydrogel material.

Materials: - Silicone Hydrogel Filcon V3 Dk 60 x 10 ⁻¹¹ 74% water content (3 months) - Hydrogel Menicon Soft 72 Dk 34 x 10 ⁻¹¹ 72% water content (6 months)	Ø _T (mm)	C _t (mm)	BC (mm)	F'v per 0.25 (D)	Edge Lift (mm)
ROSE K2 Soft™ 1 lens Daily wear 2-lens pack is available for Silicone Hydrogel	14.30 to 15.30 per 0.10 (standard: 14.80)	0.25 to 0.60 per 0.01 (standard: 0.35)	7.40 to 9.00 per 0.20	-30.00 to +30.00 Cylinders: -0.25 D to -10.00 D per 0.25 D Axis: 0° to 180° per 1°	standard 0.0 flat +1.0 steep -1.0 double flat +2.0 double steep -2.0



Menicon recommends to use MeniCare Soft or SOLOCARE AQUA® multi-purpose solution

SOFT CONTACT LENSES



VINTAGE

CONTACTLENZEN

Replacement

6
Months

Made to order lenses Material: hydrogel Mipafilcon A 72% water content - Dk 34 x 10 ⁻¹¹		Ø _T (mm) per 0,10mm	BC (mm) per 0,10mm	Power (D)	Cylinder (D) Axis (°)
Vintage Daily wear		13,50 to 15,50	7,00 to 9,50	-25,00 D to +25,00 D per 0.25 D	-
Vintage Toric Daily wear		13,50 to 15,50	7,00 to 9,50	-20,00 D to +20,00 D per 0.25 D	-0,75 to -3,00 D per 0,25 D Axis: 0° to 180° per 1°
Vintage Multi Daily wear		13,50 to 15,50	7,00 to 9,50	-25,00 D to +25,00 D per 0.25 D Addition +1,00 to +3,00 per 0,50 D	-
Vintage Multi Toric Daily wear		13,50 to 15,50	7,00 to 9,50	-20,00 D to +20,00 D per 0.25 D Addition +1,00 to +3,00 per 0,50 D	-0,75 to -3,00 D per 0,25 D Axis: 0° to 180° per 1°

Menicon recommends to use MeniCare Soft or SOLOCARE AQUA®, multi-purpose solution

Attention: Menicon recommends to replace Menicon soft 72 Toric by Vintage Toric

Replacement

1
Year

Annual lenses

Conventional lenses Material: Hydrogel Mipafilcon A 72% water content - Dk 34 x 10 ⁻¹¹		Ø _T (mm)	BC (mm) per 0,30mm	Power (D) per 0.25 D
Menicon soft 72 Daily wear		13.50	8.10 8.70	-25.00 to +25.00
		14.00	8.40 to 9.30	

Menicon recommends to use MeniCare Soft or SOLOCARE AQUA® multipurpose solution.

Attention: Menicon recommends to replace Menicon soft 72 Toric by Vintage Toric (see above this page)

SPHERIC CONTACT LENSES

Replacement 1 Year

Z MATERIAL	Z Material: Tisofcon A Dk 163 x 10 ⁻¹¹ (ISO)	Ø _T (mm)	BC (mm)	Power (D) per 0.25 (D)
	Menicon Z Daily wear or extended up to 30 days (29 nights max) Aspheric	9.00 9.20 9.40 9.60 9.80	7.20 to 8.40 per 0.10 mm	-25.00 to +25.00
	Menicon Z-α Daily wear or extended up to 30 days (29 nights max) Sphero-aspheric	8.00 to 10.00 (0.20) 10.00 10.50 11.00	6.50 to 9.00 per 0.05 mm 6.70 to 9.00 per 0.05 mm 7.30 to 9.00 per 0.05 mm 7.50 to 9.00 per 0.05 mm	-25.00 to +25.00
	Menicon EX-Z Daily wear or extended up to 30 days (29 nights max) Sphero-aspheric	8.80 9.20 9.60 10.00 10.50 ⁽¹⁾ 11.00	6.50 to 9.00 per 0.05mm 6.90 to 9.00 per 0.05mm 7.20 to 9.00 per 0.05mm	-25.00 to +25.00
	Menicon Z Comfort Daily wear with tear-reservoir zone	10.20 10.60 11.00	6.50 to 9.20 per 0.05mm	-25.00 to +25.00
	Menicon Z OMNI Daily wear Sphero-aspheric	8.80, 9.10, 9.60, 10.10, 10.50, 11.00	6.00 to 9.95 per 0.05mm	-25.00 to +25.00

EX MATERIAL	EX Material: Tolofcon A - Dk 64 x 10 ⁻¹¹ (ISO)	Ø _T (mm)	BC (mm) per 0.05 mm		Power (D) per 0.25 (D)	
			Standard	Customized manufacturing	Standard	Customized manufacturing
	Menicon EX Daily wear Sphero-aspheric	7.50 to 8.70 per 0.10 8.80 to 10.00 per 0.40 ⁽¹⁾ 10.50, 11.00 ⁽¹⁾ 11.10 to 12.00 per 0.10	6.50 to 9.00 per 0.05mm	6.00 to 6.45 9.05 to 9.40	-25.00 to +25.00	-25.25 to -30.00 +25.25 to +30.00

(1) Lenticular geometry EX 11.00mm Ø₀: 7.00mm

TORIC CONTACT LENSES

Replacement 1 Year

Z, EX MATERIAL ⁽²⁾		Ø _T (mm)	BC (mm) per 0.05mm	Cylinder (D)	Power (D) per 0.25 (D)
	Menicon Z BTC Daily wear (Back toric compensation)	9.10 9.60 10.10	7.00 to 9.95	Back cylinder (mm) : 0.20 to 1.50 per 0.05 ⁽³⁾	-25.00 to +25.00
	Menicon Z Comfort BTC Daily wear (Back toric with large diameter and tear-reservoir zone)	10.20 10.60 11.00	6.50 to 9.20	Back cylinder (mm) : 0.20 to 1.00 per 0.05 ⁽³⁾	-25.00 to +25.00
	Menicon Z BT Daily wear (Bitoric)	9.10 9.60 10.10	7.00 to 9.95	Back cylinder (mm) : 0.20 to 1.50 per 0.05 ⁽⁴⁾ Tore total (D) ≤ 10 D ⁽³⁾ Axis: 0° to 180° per 1°	-25.00 to +25.00
	Menicon Z Comfort BT Daily wear (Bitoric with large diameter and tear-reservoir zone)	10.20 10.60 11.00	7.00 to 9.95	Back cylinder (mm) : 0.20 to 1.00 per 0.05 ⁽³⁾ Tore total (D) ≤ 10 D ⁽⁴⁾ Axis: 0° to 180° per 1°	-25.00 to +25.00
	Menicon Z FT Daily wear (Front Toric)	9.10 9.60 10.10	6.50 to 9.95	Front Cylinder (D): -0.75 à -2.00 par 0.25 (D) Axis: 0° to 180° per 1° Prism ballast: 1.00 to 2.50 per 0.50 Δ (standard: 1.50)	-25.00 to +25.00

(2) Available in Dk 58 material upon request. (3) Depending on the value of the flat BC. (4) Depending on the value of the back cylinder.

Note: ØT 8.80, 10.50 and 11.00mm are available for the lenses Z BTC, Z BT and Z FT.

Menicon recommends to use MENICARE PURE or MENICARE PLUS, multi-purpose solution for daily care + PROGENT once a week.

MULTIFOCAL CONTACT LENSES

Replacement

1
Year

Z, EX MATERIAL ⁽¹⁾	SPHERICAL	SIMULTANEOUS VISION Z, EX Material or Dk 58	Ø _T (mm)	BC (mm) per 0.05mm	Tore	Power (D) per 0.25 (D)	Addition
		Menicon Z Progressive <i>Daily wear</i> (Center Distance)	9.10 9.60 10.10	6.00 to 9.95		-25.00 to +25.00	<i>Addition: +1.00 to +3.00 per 0.50 D</i>
		Z Comfort Progressive <i>Daily wear</i> (Center Distance, large diameter, tear-reservoir zone)	10.20 10.60 11.00	6.50 to 9.20		-25.00 to +25.00	<i>Addition: +1.00 to +3.00 per 0.50 D</i>
		Z Comfort Progressive Near <i>Daily wear</i>	10.20 10.60 11.00				
		Z Progressive Near <i>Daily wear</i> (Center Near)	9.10 9.60 10.10	6.00 to 9.95		-25.00 to +25.00	<i>Addition: +1.00 to +3.00 per 0.50 D</i>
		Z Progressive + <i>Daily wear</i> (Accelerated near vision progression)	9.10 9.60 10.10				
	BACK TORIC	Z Progressive BTC <i>Daily wear</i>	9.10 9.60 10.10	7.00 to 9.95	<i>Back cylinder (mm): 0.20 to 1.50 per 0.05⁽²⁾</i>	-25.00 to +25.00	<i>Addition: +1.00 to +3.00 per 0.50 D</i>
		Z Comfort Progressive BTC <i>Daily wear</i>	10.20 10.60 11.00	6.50 to 9.20	<i>Back cylinder (mm): 0.20 to 1.00 per 0.05⁽²⁾</i>		
		Z Progressive + BTC <i>Daily wear</i>	9.10 9.60 10.10	7.00 to 9.95	<i>Back cylinder (mm): 0.20 to 1.50 per 0.05⁽²⁾</i>		
	BITORIC	Z Progressive BT <i>Daily wear</i>	9.10 9.60 10.10	7.00 to 9.95	<i>Back cylinder (mm) : 0.20 to 1.50 per 0.05 Tore total (D) ≤ 10 D⁽³⁾ Axis: 0° to 180° per 1°</i>	-25.00 to +25.00	<i>Addition: +1.00 to +3.00 per 0.50 D</i>
		Z Comfort Progressive BT <i>Daily wear</i>	10.20 10.60 11.00	7.00 to 9.95	<i>Back cylinder (mm) : 0.20 to 1.00 per 0.05 Tore total (D) ≤ 10 D⁽³⁾ Axis: 0° to 180° per 1°</i>		
		Z Progressive + BT <i>Daily wear</i>	9.10 9.60 10.10	7.00 to 9.95	<i>Back cylinder (mm) : 0.20 to 1.50 per 0.05 Tore total (D) ≤ 10 D⁽³⁾ Axis: 0° to 180° per 1°</i>		
	ALTERNATING VISION		Ø _T (mm)	BC (mm)		Power (D) per 0.25 (D)	
	Menicon Z Executive <i>Daily wear</i> (Segmented)		9.10 9.60 10.10	6.50 to 9.00 per 0.05mm <i>Prism ballast: 1.00 to 2.50 per 0.50 Δ (standard: 2.00). Position of segment: +0,5 mm / -0,5 mm</i>		-25.00 to +8.00 <i>Addition: +1.00 to +3.00 per 0.50 D</i>	
	Menifocal Z <i>Daily wear or extended up to 30 days (29 nights max)</i> (Concentric center distance) (Z material only)		9.60 9.80	7.00 to 9.00 per 0.10mm		-13.00 to +5.00 <i>Addition: +1.00 to +2.50 per 0.50 D</i>	

(1) Available in Dk 58 material upon request. (2) Depending on the value of the flat BC. (3) Depending on the back cylinder.

Note: ØT 8.80, 10.50 and 11.00mm are available for Z Progressive, Z Progressive Near, Z Progressive +, Z Progressive BTC, Z Progressive + BTC, Z Progressive BT, Z Progressive + BT, Z Executive.


Why should you recommend to wear contact lenses in Z material?

From children to presbyopia age, Z material can be the first choice, granting not only excellent visual quality but also daily comfort.

Thanks to its high DK (163×10^{-11}), the lenses allow the cornea to breathe through the material, and Menicon designs enable the good tears exchange.

Ideal for someone who has sensitive eyes, with dry eye syndrome, allergies or with intolerance to soft lenses.

We can start anytime: children, even babies wear rigid gas permeable lenses. Menicon offers wide range of designs with Z material (or EX) which enables adaptation to different profiles and needs. Wearing the lenses with Z material to have good near vision as soon as the first sign of presbyopia appears, or to correct the corneal astigmatism and certain forms of irregular corneas... Several solutions can be proposed by customizing to the needs of visual acuity of each person (according to its profession, activities, life-style, etc...).





ORTHOKERATOLOGY

Replacement

1 Year

Z MATERIAL	Corneal reshaping overnight lens. Dk: 163×10^{-11} (ISO)	$\varnothing_T$ (mm)	BC (mm)	Indications
	Menicon Z Night Overnight wear		7.00 to 10.00 mm per 0.05 mm <i>Tangential angle: 46° to 63° (in 1° steps)</i> <i>Sagittal depth: 0.30 to 1.00 mm (in 0.01 mm steps)</i>	Myopia ≤ -4.00 D Direct corneal astigmatism ≤ 1.50 D Inverse corneal astigmatism ≤ 0.75 D
	Menicon Z Night Toric Overnight wear	10.20 10.60 11.00	7.00 to 10.00 mm per 0.05 mm <i>Tangential angle: 46° to 63° (in 1° steps)</i> <i>In the flat meridian (FL) and in the steep meridian (ST)</i> <i>Sagittal depth: 0.30 to 1.00mm (in 0.01 mm steps)</i> <i>In the flat meridian (FL) and in the steep meridian (ST)</i>	Myopia ≤ -4.00 D Direct corneal astigmatism ≤ -2.50 D Inverse corneal astigmatism ≤ -1.50 D

Menicon recommends to use MENICARE PURE or MENICARE PLUS multi-purpose solution for daily care + PROGENT once a week.

SCLERAL LENS



TIMEXL RANGE

Replacement

1 Year

Z MATERIAL	Scleral lens with bitangential periphery in Z material. Daily wear	$\varnothing_T$ (mm)	Tangent angle	Sagittal height	BC (mm)	Power per 0.25 (D)
	TimeXL BT* * Spherical design upon request, without front cylinder	14.00 to 17.50 per 0.05mm	28° to 50° per 2°	3000 to 4600 μ m	7.40 to 9.20 per 0.20mm	-25.00 to +25.00 D per 0.25 D Cylinders: -0.75 D to -3.00 D per 0.25 D Axis: 0° to 180° per 1°

Menicon recommends to use MENICARE PURE or MENICARE PLUS multi-purpose solution for daily care + PROGENT once a week. Use special lens case for large diameter lenses. Trial lens set is available

STARTER KIT LARGE DIAMETER RGP LENSES

New

- 1 MENICARE PURE 70ml
- + 2 doses A & B PROGENT
- + 1 lens case large diameter PROGENT
- + 1 lens case large diameter
- + 1 suction cup DMV® MAXI CUP™ (to put)
- + 1 suction cup DMV® Ultima™ (to remove)

Order unit: 5 kits

Product code: 2321





ROSE K LENSES

Replacement

1
Year

Z, EX MATERIAL OR Dk 65		Ø _T (mm)	BC (mm)	Power per 0.25 (D)	Edge Lift (mm)
	ROSE K2 Daily wear (Aberration control) ROSE K Daily wear	7.50 to 11.00 per 0.10	4.30 to 8.80 per 0.05 <i>*Available variations: +3.0 to -1.3 per steps 0.1</i>	+ 40.00 to -40.00	standard 0.0 flat +1.0 steep -0.5 others*
	ROSE K2 POST GRAFT Daily wear	9.00 to 12.50 per 0.10 (standard: 10.40)	5.20 to 11.60 per 0.05 <i>Available variations: +3.0 to -3.0 per steps 0.5</i>	+ 40.00 to -40.00	standard 0.0 flat +0.6 - steep -0.5 double flat +1.2 double steep -1.0
	ROSE K2 IC (Irregular Cornea) Daily wear	9.00 to 12.50 per 0.10 (standard: 11.20) 11.20, r0 ≥ 6.00 12.00, r0 ≥ 6.70 <i>Available variations: +3.0 to -3.0 per steps 0.5</i>	5.30 to 10.00 per 0.05	+ 40.00 to -40.00	standard 0.0 flat +0.6 - steep -0.5 double flat +1.2 double steep -1.0
	ROSE K2 NC (Nipple Cone) Daily wear	7.60 to 10.00 per 0.10	4.00 to 8.10 per 0.05 <i>Available variations: +4.0 to -1.5 per steps 0.1</i>	+15.00 to -40.00	standard 0.0 flat +1.0 steep -0.5
Z MATERIAL	ROSE K2 XL (Z material only) Daily wear CLEAR tint is available	13.00 to 16.60 per 0.10 (standard: 14.60)	5.60 to 9.60 per 0.05 <i>Available variations: +6.0 to -4.0 per steps 0.5</i>	+40.00 to -40.00	standard 0.0 flat +1.0 steep -1.0 double flat +2.0 double steep -2

ROSE K2, ROSE K2 Post Graft, ROSE K2 IC

ROSE K2 XL Daily wear

Option toric or asymmetric (CE 1639)

- Toric peripheral curves: 0.40 to 2.00mm.
- Asymetric corneal technology (ACT): 0.40 to 1.50mm (in 0.10mm steps). (0.80 to 2.00mm for ROSE K2 XL).
- Toric back, front and bi-toric surfaces, axis: 0° to 180° (in 1° steps).

Menicon recommends to use MENICARE PURE or MENICARE PLUS multi-purpose solution for daily care + PROGENT once a week.

For material availability, other possibilities of parameters and asymmetric design, please contact our Technical Assistance:
atechnique@menicon.fr - ROSE K2 XL: Z Material only. Trial lens set is available for each range: K2, Post Graft, IC, NC or XL.

SOFT LENS FOR IRREGULAR CORNEAS

ROSE K2 Soft is a daily wear soft lens for irregular corneas, 3 months replacement lens manufactured from Silicon Hydrogel material, or 6 months or annual replacement lens manufactured from Hydrogel material (page 6).

RIGID GAS PERMEABLE LENSES - OTHER DESIGNS

Z MATERIAL OR EX		Ø _T (mm)	BC (mm) per 0.05	Power per 0.25 (D)
	Menicon Plateau Daily wear	8.00 to 11.00 (standard: 10.00) per 0.10	7.03 to 10.23 (48.00 to 33.00 D) r ₁ 1.00 to 7.00 D per 1.00 D	+20.00 to -20.00
	Menicon S9 Daily wear	8.50 8.80 9.00 9.20 9.50	6.50 to 9.00	+25.00 to -25.00
	Menicon KRC 1 to 12 Daily wear	7.90 8.00 8.10	5.15 to 6.90	-0.25 to -25.00
	Menicon K3 Daily wear	8.50	5.00 to 6.90	-0.25 to -25.00
	Menicon TCC6 Daily wear	9.20 9.30 9.50 10.00 10.50	6.50 to 9.00	+25.00 to -25.00
	Menicon B4P Daily wear	10.00 10.50 11.00 11.50 12.00	6.50 to 9.00 7.20 to 9.00 7.25 to 9.00	+25.00 to -25.00
	Menicon B4PM Daily wear	11.00 11.50 12.00	7.20 to 9.00 7.25 to 9.00	+25.00 to -25.00
	Menicon MX Daily wear	8.50 8.80 9.00	6.50 to 9.00	+25.00 to -25.00
Other designs, please contact our Technical Assistance				

Menicon recommends to use MENICARE PURE or MENICARE PLUS multi-purpose solution for daily care + PROGENT once a week.

LENS CARE SOLUTIONS

FOR ALL RIGID GAS PERMEABLE CONTACT LENSES

MENICARE PURE 250ML

Multipurpose solution (daily use)
250ml
+ lens case

Order unit: 10

Product code: 2312



PROGENT (5 treatments)

Intensive cleaner for protein
and for disinfection
(weekly use)

Box of 5 doses A
+ 5 doses B
+ 1 PROGENT vial

Order unit: 10

Product code: 1156



MENICARE PURE & PROGENT

STARTER KIT

MENICARE PURE 70ml
1 bottle
+ 1 lens case, 1 dose of
PROGENT A&B
+ 1 PROGENT vial

Order unit: 10

Product code: 2316



MENICARE PURE 70ML

Multipurpose solution (daily use)
70ml
+ lens case

Order unit: 10

Product code: 2315



MENICARE PLUS 250ML

Multipurpose solution
(daily use)
250ml
+ lens case
(Also available in Starter Kit :
Menicare Plus 50ml
+ 1 lens case
+ 1 dose of Progent A & B
+ 1 Progent lens vial)

Order unit: 10

Product code: 2112



SPRAY & CLEAN (CLEANER)

Preservative free solution
15ml
(cleaning before wear)
Antilipid action
(daily use)

Order unit: 10

Product code: 1254



MENILAB (trial lenses)

For professional use
Disinfection
of RGP trial lenses
250ml

Order unit: 1

Product code: 1852



ACCESSORIES FOR ALL RIGID GAS PERMEABLE CONTACT LENSES

RGP lens case (Menicase)

individual lens case

Order unit: 10

Product code : 5400 C



Suction cup (Menicup)

Order unit: 10 Product
code : 5500 C



Travel Kit (Menibox)

Case with mirror,
2 lens cases, 1 suction cup

Order unit: 1

Product code : 5800 C



FOR SEMI-SCLERAL CONTACT LENSES (ROSE K2 XL, TIME XL)

Suction DMV MaxiCup™

(Insertion)

Order unit: 3

Product code : DMV CUP



Suction cup DMV® Ultima™

(Removal)

Order unit: 10

Product code : DMV ULTIMA



Large diameter RGP lens case (PROGENT treatment)

Order unit: 1

Product code : ETUI LRPO



Flat lens case

(For large diameter lenses, storage case)

Order unit: 1

Product code : ETUI GDDIA



LENS CARE SOLUTIONS FOR ALL SOFT CONTACT LENSES

MENI CARE SOFT

Available
from
April 2020



New family:

- New design
- New bottles of 250ml and 70ml



■ MENICARE SOFT 250ML

Multipurpose solution (daily use)
250ml
+ lens case

Order unit: 10

Product code: 2208



New cap*



New opening system
(screw the cap in clockwise direction)*

■ MENICARE SOFT 70ML

Multipurpose solution (daily use)
70ml
+ lens case

Order unit: 10

Product code: 2209



■ TRAVEL KIT MENICARE SOFT 3x70ML

Multipurpose solution
(daily use)
3 x 70ml
+ 3 lens cases

Order unit: 10

Product code: 2211



■ MENICARE SOFT 360ML

Multipurpose solution (daily use)
360ml
+ lens case

Order unit: 10

Product code: 2202

MENICARE SOFT in 360ml new design will be available during the year 2020



*Only for 70ml bottle and 250ml bottle

LENS CARE SOLUTIONS FOR ALL SOFT CONTACT LENSES

SOLOCARE AQUA®
by  **Menicon**



Multipurpose solution
with MicroBlock® lens case

■ SOLOCARE AQUA® 360ML

Multipurpose solution (daily use)
360ml
+ lens case

Order unit: 10

Product code: 6162



Only available in EEC

■ TRAVEL KIT SOLOCARE AQUA® 90ML

Multipurpose solution
(daily use)
90ml
+ lens case
+ mirror

Order unit: 12

Product code: 6207



Only available in EEC

■ SOLOCARE AQUA® 90ML

Multipurpose solution (daily use)
90ml
+ lens case

Order unit: 12

Product code: 6158



Only available in EEC

ACCESSORIES FOR SOFT CONTACT LENSES

■ HOMECASE

Storage accessory
30 "Flat Pack" (R) / 30 "Flat Pack" (L)

Order unit: 1

Product code: HOME CASE



■ TWEEZERS (Meniholder)

Order unit: 10

Product code: 5900 C





Fit Gas Permeable & Soft Contact Lenses
quickly & easily


Free download:
<http://www.menicon.fr/pro/fitting-software/easyfit/>

WORKS WITH OR WITHOUT TOPOGRAPH

Four functions:

- 1 Recommended suitable Menicon lenses based on Topography data
 - RGP spherical, toric and multifocal
 - TimeXL Range
 - ROSE K2 Range
 - Menicon Z Night, Menicon Z Night Toric
 - Miru 1day Flat Pack, Miru 1day UpSide
 - Menicon PremiO Range
 - Miru 1month Menicon Range
 - Individual Range
 - Vintage Range

- 2 Order directly from Easyfit
- 3 Get technical advice more easily
 - Send the patient data file just with one click
- 4 Automatic sending of orders by mail



Realistic and dynamic fluorescein
images

With Menicon design lenses

Compatible with more than 27 topographers

BON Eyetop	Phoenix
Cassini	Rodenstock DNEye Scanner
EyeTop 2005	Shin-Nippon CT-1000
Haag-Streit	Sirius
Keratron	Tomey RT-7000
Essilor Anaeyes	Tomey TMS-2 / TMS-3 / TMS-4 / TMS-5
Medmont (version 3)	Topcon 7000
Medmont (version 4)	Topcon 8100 (CA-100 software)
Medmont (version 5)	Topcon 8100 (Color Mapping software)
Medmont (version 6)	Topcon CA-100 / CA-200 / CA-800
Nidek OPD-Scan III	Visionix VX-120
Oculus Easygraph	Zeiss i-Profilor
Oculus Keratograph	Zeiss i-Profilor Plus
Oculus Pentacam	

Access to technical documentations

Information: www.e-menicon.fr

Keratron Piccolo

Topographer and Aberrometer

Accessories upon request

Extremely reliable and precise
Online support



MENICON GLOBAL EXPERTISE TO FIT

Even the most complex cases

Simply find the solution
From the wide range of
contact lenses

EASYFIT SOFTWARE

To download Easyfit :
<http://www.menicon.fr/pro/fitting-software/easyfit/>
or contact customer service : serviceclients@menicon.fr

TECHNICAL ASSISTANCE

Monday to Thursday from 9:00 am to 6:00 pm
Friday from 9:00 am to 5:00 pm
Phone: +33 1 85 73 13 12
atechnique@menicon.fr

TRAINING PROGRAM

Contact us for more information
campus@menicon.fr

CUSTOMER SERVICE

Monday to Thursday from 9:00 am to 6:00 pm
Friday from 9:00 am to 5:00 pm
Phone: +33 1 85 73 13 11 - Fax: +33 1 85 73 13 97
serviceclients@menicon.fr

TRIAL LENSES, STARTER KITS FOR LENS CARE

Contact customer service:
serviceclients@menicon.fr

