



# MENICON CO., Ltd., Japan **SERVIS**

## MEDECO & AREDOVALE WEB STRÁNKY

• www.menicon.sk • www.menicon.cz

## ZÁKAZNÍCKY SERVIS, OBJEDNÁVKY

### Slovenská Republika

MEDECO spol. s r.o. M: +421 948 771 277  
Bezekova 2034/ 21 M: +421 905 946 919  
Bratislava 841 02 Z: +421 800 328 800  
hummel@medeco.sk sino@medeco.sk info@medeco.sk

### Česká Republika

AREDOVALE s.r.o.  
Karla Engliše 3201/ 6  
Praha 5 - Smíchov, 150 00  
M: +421 905 946 919 +420 800 722 227  
info@opto-shop.cz aredovale@aredovale.cz

*Od chvíle, kdy otevřeme ráno oči, dokud je v noci znovu nezavřeme, naše dny a náš život - jsou určeny tím, co vidíme. Proto jsme si vybrali Miru - japonské slovo, které znamená "vidět" - jako globální značku naší nové vizionářské řady kontaktních čoček.*

*Spolu s dalšími produkty výrobní řady Miru se navždy mění způsob, jakým vidíme svět.*




*Od okamihu, keď ráno otvoríme oči, až kým ich opäť v noci nezatvoríme, naše dni a naše životy - sú určované tým, čo vidíme. Preto sme si vybrali Miru - japonské slovo, ktoré znamená "vidieť" - ako globálnu značku nášho nového vizionárskeho radu kontaktných šošoviek.*



*Spolu s ostatnými produktmi výrobnej rady Miru, sú stvorené navždy zmeniť spôsob nášho videnia sveta.*

## MÄKKÉ (SOFT) KONTAKTNÉ ŠOŠOVKY/MĚKKÉ ( SOFT ) KONTAKTNÍ ČOČKY

Doba výmeny	1 Day		$\varnothing_T$ (mm)	$r_0$ (mm)	Power (D)	Cylinder (D) Axis (°)
Miru 1day Menicon Flat Pack Box of 30 lenses Material: Hioxifilcon A (HEMA - GMA) 57 % water content Dk/e : 25x10 <sup>-9</sup>			14.20	8.60	+4.00 to -6.00 D in 0.25 D steps -6.50 to -10.00 D in 0.50 D steps	
Miru 1day UpSide Box of 30 lenses Material: Midafilcon A 59 % water content Dk/ t: 91 x10 <sup>-9</sup> Central thickness 0,07 mm -3,00 D			14.20	8.40	-0.25 to -6.00 D in 0,25 D steps -6.00 to -10.00 D in 0.50 D steps	

# MÄKKÉ (SOFT) KONTAKTNÉ ŠOŠOVKY/MĚKKÉ ( SOFT ) KONTAKTNÍ ČOČKY

Doba výmeny <b>1</b> Mesiac		$\varnothing_T$ (mm)	$r_0$ (mm)	Power (D)	Cylinder (D) Axis (°)
<b>Miru 1month Menicon</b> Box of 6 lenses Material: Silicone hydrogel Asmofilcon A 40% water content - Dk 129 x10 <sup>-11</sup> Dk/e : 161x10 <sup>-9</sup>		14.00	8.30 8.60	+6.00 to -6.00 D in 0.25 D steps  -6.50 to -13.00 D in 0.50 D steps	-
<b>Miru 1month Menicon For Astigmatism</b> Box of 6 lenses Material: Silicone hydrogel Asmofilcon A 40% water content - Dk 129 x10 <sup>-11</sup> Dk/e : 161x10 <sup>-9</sup>		14.00	8.60	00.00 to -6.00 D in 0.25 D steps  -6.50 to -10.00 D in 0.50 D steps	Cylinders -0.75 D -1.25 D -1.75 D Axis 10°- 180° v 10°steps
<b>Miru 1month Menicon Multifocal</b> Box of 6 lenses Material: Silicone hydrogel Asmofilcon A 40% water content - Dk 129 x10 <sup>-11</sup> Dk/e : 161x10 <sup>-9</sup>		14.20	8.60	+6.00 to -6.00 D in 0.25 D steps  -6.00 to -13.00 D in 0.50 D steps	Add: LOW (+1.00) HIGH (+2.00)

Doba výmeny <b>3</b> Mesiac		$T$ (mm)	$r_0$ (mm)	Power (D)	Cylinder (D) Axis (°)
<b>Range Individual</b> 1 šošovka, alebo 1 balenie 2 ks šošovky <b>Individual</b> <b>Individual Toric</b> <b>Individual Progressive 2*</b> (Central Distance) <b>Individual Progressive 2 Toric</b> <b>Individual Progressive 2 ear</b> (Central near) <b>Individual Progressive 2 ear Toric</b>		13.00 to 15.50 per 0.10 mm	7.40 to .50 per 0.10 mm	-25.00 to 25.00 D in 0.25 D steps Addition 1.00 to 3.00 in 0.50 D steps	Cylinders -0.75 D to -5.00 D -0.75 D to -3.00 D (Prog. 2 Toric and Prog. 2 Near Toric) in 0.25 D steps Axis from 0° to 180° per 1°
		Material: Silicone hydrogel Filcon V3 72% water content - Dk 60 x10 <sup>-11</sup> *Individual Progressive is also available			
<i>Keratoconus contact lenses and special cases</i>					
<b>ROSEK 2 Soft</b> 1 lens or 1 box of 2 lenses with same parameters		(Ref page)			

## SOLOCARE AQUA®



### Tri Pack SOLOCARE AQUA®

3 x 360ml + 3 lens case

Order unit: 6



### Duo Pack SOLOCARE AQUA®

2 x 360ml + 2 lens case

Order unit: 12



### SOLOCARE AQUA®

360ml + lens case

Order unit: 10



### SOLOCARE AQUA®

90ml + lens case

Order unit: 12



### Travel kit SOLOCARE AQUA®

90ml + lens case + mirror

Order unit: 12

## MENICARE SOFT



### Pack MeniCare Soft

2 x 360ml + 2 lens case

Order unit: 10



### MeniCare Soft

360ml + lens case

50ml + lens case / 70ml + lens case\*

Order unit: 10

\*the volume of small bottle will be changed from 50ml to 70ml

## KERATOCONUS IRREGULAR CORNEA

### SOFT LENS FOR IRREGULAR CORNEAS

ROSE K2 Soft is a daily wear soft lens for irregular corneas, 3 month replacement lens manufactured from Silicon Hydrogel material, or 6 month or annual replacement lens manufactured from Hydrogel material.



	$\varnothing_T$ (mm)	$e_C$ (mm)	$r_0$ (mm)	F'v in 0.25 (D) steps	Edge Lift (mm)
<b>ROSE K2 Soft</b> Materials : -Silicone hydrogel Lagado (LSiH), Dk 49 49% water content  -Menicon soft 72 Hydrogel Dk 34 72% water content	14.30 to 15.30 in 0.10 steps (standard : 14.80)	0.25 to 0.60 in 0.10 steps (standard : 0.35)	7.40 to 9.00 in 0.20 steps	-30.00 to +30.00  Cylinders -0.25 D to -10.00 D in 0.25 D steps  Axis from 0° to 180° per 1°	standard 0.0 flat +1.0 steep -1.0 double flat +2.0 double steep -2.0

## SFERICKÉ

	$\varnothing_T$ (mm)	$r_0$ (mm) in 0.05 mm steps		Power (D) in 0.25 (D) steps		
		Standard	Non-standard	Standard	Non-standard	
<b>Z MATERIAL</b>	Menicon Z Aspheric Material : Tisilfocon A Dk: 163 (ISQ)	9.00 9.20 9.40 9.60 9.80	7.20 7.30 to 7.50 7.60 to 7.80 7.90 to 8.10 8.20 to 8.40 in 0.10 mm steps	$\varnothing_T$ from 9.00 to 9.80 for 7.20 to 8.60 in 0.10 mm steps	+5.00to-10.00	+5.25to+25.00 -10.25 to -25.00
	Menicon Z□ Sphero-aspheric	8.00 to 9.00 in 0,20 mm steps 8.80 9.20 9.40 9.60 9.80 10.00 10.50 11.00	7.20 to 8.40 in 0.05 mm steps	6.50 to 9.00 in 0.05 mm steps  6.70 to 9.00 7.30 to 9.00 7.50 to 9.00	+5.00 to -10.00	+5.25to+25.00 -10.25 to -25.00
	Menicon EX-Z Sphero-aspheric	8.80 9.20 9.60 10.00 10.50 11.00*	7.20 to 8.40 in 0.05 mm steps	6.50 to 8.60 6.50 to 8.80 6.50 to 9.00 6.90 to 9.00  7.20 to 9.00	+5.00 to -10.00	+5.25to+25.00 -10.25 to -25.00
	Menicon Z Comfort with tear-reservoir zone	10.20 10.60 11.00	7.20 to 8.40 in 0.05 mm steps	6.50 to 9.00 in 0.05 mm steps	+5.00 to -10.00	+5.25to+25.00 -10.25 to -25.00
	Menicon Z OMNI Sphero-aspheric (Z or EX material)	9.10 9.60 10.10	7.20 to 8.40 mm	6.50 to 9.00 in 0.05 mm steps	+5.00 to -10.00	+5.25to+25.00 -10.25 to -25.00
<b>EX MATERIAL</b>	Menicon EX Sphero-aspheric Material: Tolofoccon A DK: 64 (ISQ)	8.80 9.20 9.60 10.00 10.50 11.00*	7.20 to 8.40 in 0.05 mm steps	6.50 to 8.60 6.50 to 8.80 6.50 to 9.00 6.90 to 9.00  7.20 to 9.00	+5.00to-10.00	+5.25to+25.00 -10.25 to -25.00

\* Lenticular geometry EX 11.00 mm  $\varnothing_0$ : 7.00 mm

## TORICKÉ

	$\varnothing_T$ (mm)	$r_0$ (mm) in 0.05 mm steps		Power (D) in 0.25 (D) steps		
		Standard	Non-standard	Standard	Non-standard	
<b>Z OR EX MATERIAL</b>	Menicon Z BTC (Back toric compensation)	9.10 9.60 10.10	7.20 to 8.40 mm	7.00 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00
	Menicon Z Comfort BTC (Back toric with large diameter and tear-reservoir zone)	10.20 10.60 11.00	7.20 to 8.40 mm	6.50 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00
	Menicon Z BT (Bitoric)	9.10 9.60 10.10	7.20 to 8.40 mm	7.00 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00
	Menicon Z Comfort BT (Bitoric with large diameter and tear-reservoir zone)	10.20 10.60 11.00	7.20 to 8.40 mm	7.00 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00
	Menicon Z FT (Front Toric)	9.10 9.60 10.10	7.20 to 8.40 mm	6.50 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00

\*Other parameters might be available upon request, please contact our Technical Assistance

MULTIFOKÁLNE

Z OR EX MATERIAL	Simultaneous vision	$\varnothing_T$ (mm)	$r_0$ (mm) per 0.05 mm		Power (D) per 0.25 (D)	
			Standard	Non-standard	Standard	Non-standard
SPHERIC	Menicon Z Progressive (Center Distance)	9.10	7.20 to 8.40 mm	6.50 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00
		9.60				
		10.10				
	Z Comfort Progressive (Center Distance, large diameter, tear-reservoir zone)	10.20	7.20 to 8.40 mm	6.50 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00
		10.60				
		11.00				
Z Comfort Progressive Near	10.20	Addition : +1.00 to +3.00 in 0.50 D steps				
	10.60					
	11.00					
Z Progressive Near (Center Near)	9.10	7.20 to 8.40 mm	6.50 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00	
	9.60					
	10.10					
Z Progressive $\phi$ (accelerated near vision progression)	9.10	Addition : +1.00 to +3.00 in 0.50 D steps				
	9.60					
	10.10					
BACK TORIC	Z Progressive BTC	9.10	7.20 to 8.40 mm	7.00 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00
		9.60				
		10.10				
	Z Comfort Progressive BTC	10.20	7.20 to 8.40 mm	6.50 to 9.00 mm	+5.00 to -10.00	
		10.60				
		11.00				
Z Progressive + BTC	9.10	7.20 to 8.40 mm	7.00 to 9.00 mm	+5.00 to -10.00		
	9.60					
	10.10					
BITORIC	Z Progressive BT	9.10	7.20 to 8.40 mm	7.00 to 9.00 mm	+5.00 to -10.00	+5.25 to +25.00 -10.25 to -25.00
		9.60				
		10.10				
	Z Comfort Progressive BT	10.20	7.20 to 8.40 mm	7.00 to 9.00 mm	+5.00 to -10.00	
		10.60				
		11.00				
Z Progressive + BT	9.10	Addition : +1.00 to +3.00 in 0.50 D steps Back cylinder (mm) : 0.20 to 1.50 in 0.05 steps (except for Z Comfort Prog. BT : 0.20 to 1.00) Front cylinder (D) : -0.50 to -2.00 in 0.25 steps Axis: 0° to 180° in 1° steps				
	9.60					
	10.10					
Alternating vision		$\varnothing_T$ (mm)	$r_0$ (mm) per 0.05 mm		Power (D) per 0.25 (D)	
			Standard	Non-standard	Standard	Non-standard
Menifocal Z (Z material only)	9.60	7.20 to 8.40 mm per 0.05 mm	7.00 to 9.00 mm	+5.00 to -10.00	-10.25 to -13.00	
	9.80					
Menicon Z Executive (segmented)	9.10	7.20 to 8.40 mm per 0.05 mm	6.50 to 9.00 mm	+5.00 to -10.00	-10.25 to -25.00	
	9.60					
	10.10					

Upon request, lenses are available with Dk 58 material (except for Menifocal Z)

\*Other parameters might be available upon request, please contact our Technical Assistance

ORTOKERATO

Z MATERIAL	$\varnothing_T$ (mm)	$r_0$ (mm) 7.00 to 10.00 mm per 0,05 mm
Menicon Z Night		Tangential angle: 46° to 63° (in 1° steps) Sagittal depth: 0.30 to 1.00 mm (in 0.01mm steps)
	10.20	
Menicon Z Night Toric	10.60	Tangential angle: 46° to 63° (in 1° steps) In the flat meridian (FL) and in the steep meridian (ST)
	11.00	Sagittal depth: 0.30 to 1.00mm (in 0.01mm steps) In the flat meridian (FL) and in the steep meridian (ST)

# SOFTWARE

## FOR EYE CARE PROFESSIONALS

**easyfit™**

Fit Gas Permeable & Soft Contact Lenses  
quickly & easily

To get more information: Contact your Menicon regional sales manager

### Three functions:

- 1 Recommend suitable lenses based on Topography data
  - RGP spherical, toric, multifocal
  - ROSE K2 Range
  - Menicon ZNight, ZNight Toric
  - Miru 1day Menicon Flat Pack
  - Menicon PremiO, PremiO Toric
  - Miru 1month Menicon
  - Miru 1month Menicon For Astigmatism
  - Individual Range
- 2 Order directly from easyfit
- 3 Get technical advice more easily
  - Send the patient data file just with one click

### Compatible with more than 25 topographers

BON Eyetop  
Cassini  
EyeTop 2005 Haag-Streit  
Keratron  
Essilor Anaeyes  
Medmont (version 3)  
Medmont (version 4)  
Medmont (version 5)  
Medmont (version 6)  
Nidek OPD-Scan III  
Oculus Easygraph  
Oculus Keratograph  
Oculus Pentacam

Phoenix  
Rodestock DNEye Scanner  
Shin-Nippon CT-1000  
Tomey RT-7000  
Tomey TMS-2 / TMS-3 / TMS-4  
Topcon 7000  
Topcon 8100 (CA-100 software)  
Topcon 8100 (Color Mapping software)  
Topcon CA-100 / CA-200 / CA-800  
Visionix VX-120  
Zeiss i-Profiler  
Zeiss i-Profiler Plus

## STAROSTLIVOST' O TVRDÉ KONTAKTNÉ ŠOŠOVKY / PÉČE O PEVNÉ KONTAKTNÍ ČOČKY

### MENICARE PURE (multi-purpose solution)

MeniCare Pure  
Multi-purpose solution  
(daily use)  
250ml + lens case  
70 ml+ lens case  
Order unit: 10



### PROGENT (protein remove)

New packaging  
Menicon Progent (5 treatment)  
Intensive cleaner for protein and for disinfection  
(weekly use)  
Box of 5 doses A + 5 dose B + 1 Progent vial  
Order unit: 10



### NEW INSEPARABLE PACKS (MENICARE PURE & PROGENT)

Starter Kit  
MeniCare Pure & Progent  
MeniCare Pure 70ml 1 bottle  
+ 1 lens case, 1 dose of Progent A&B  
+ 1 Progent vial  
Order unit: 10



# KERATOCONUS IRREGULAR CORNEA

## ROSE K LENSES



Z OR EX MATERIAL		$\varnothing_T$ (mm)	$r_0$ (mm)	F'v per 0.25 (D)	Edge Lift (mm)
	ROSE K2 (Aberration contro) ROSE K	7.90 to 10.40 in 0.10 steps (standard : 8.70)	4.30 to 8.55 in 0.05 steps	+ 25.00 to -25.00 <sup>1</sup> from +3 decreased to -1.3 increased in 0.10 increments	standard 0.0 flat +1.0 steep -0.5 others <sup>1</sup>
	ROSE K2 POST GRAFT	9.40 to 12.00 in 0.10 steps (standard : 10.40)	5.70 to 9.30 in 0.10 steps	+25.00 to -25.00	standard 0.0 flat +0.6 steep -0.5 double flat +1.2 double steep -1.0
	ROSE K2 IC	9.40 to 12.00 in 0.10 steps (standard : 11.20 11.20, $r_0$ $\square$ 6.00 12.00, $r_0$ $\square$ 6.70)	5.70 to 9.30 in 0.10 steps	+25.00 to -25.00	standard 0.0 flat +0.6 steep -0.5 double flat +1.2 double steep -1.0
	ROSE K2 NC	7.60 to 9.00 in 0.10 steps (standard : 8.30)	4.30 to 7.60 in 0.05 steps	+10.00 to -35.00 in 0.25 steps <sup>2</sup> Available variations -1.5 to +3.00 in 0.1 steps	standard 0.0 steep -0.5 flat+1.0 others <sup>2</sup>
	ROSE K2 XL (Z material only)	13.60 to 15.60 in 0.10 steps (standard : 14.60)	5.80 to 8.40 in 0.10 steps	+15.00 to -30.00 in 0,25 steps <sup>3</sup> from +3.00 decreased to -3.00 increased in 0.50 increments	standard 0.0 flat+1.00 steep -1.00 double flat+2.0 double steep -2,00 others <sup>3</sup>

ROSE K2, ROSE K2 Post Graft, ROSE K2 IC

ROSE K2 XL

) Option toric or asymmetric:

- Toric peripheral curves: 0.40 to 2.00mm
- Asymetric corneal technology (ACT): 0.40 to 1.50mm (in 0.10mm steps)  
(0.80 to 2.00mm for ROSE K2 XL)
- Toric back, front and bi-toric surfaces, axis : 0° to 180° (in 1° steps)

## CENNÍK PRODUKTOV

ceny bez DPH

Miru 1day UpSide 30 ks	16,99 EUR	479 Kč
Miru 1day Menicon Flat Pack	13,50 EUR	359 Kč
Miru 1month Menicon	17,90 EUR	469 Kč
Miru 1month Menicon for Astigmatism	21,90 EUR	579 Kč
Miru 1month Menicon Multifocal	32,90 EUR	819 Kč
SOLOCARE AQUA 360 ml	8,00 EUR	195 Kč
SOLOCARE AQUA 90 ml	3,04 EUR	72 Kč
SOLOCARE AQUA 2x360 ml	14,45 EUR	350 Kč
Indivisual Spheric 1 ks	44,00 EUR	1.120 Kč
Indivisual Toric 1 ks	54,00 EUR	1.400 Kč
Indivisual Progressive 2 1 ks	48,00 EUR	1.220 Kč
MeniCare Soft 360 ml	6,90 EUR	179 Kč
MeniCare Soft 50 ml	2,90 EUR	75 Kč
RGP sférické 1 ks	50,00 EUR	1.280 Kč
RGP tórické 1 ks	90,00 EUR	2.295 Kč
ROSE K2/ K2 XL	100,00 / 150,00 EUR	2.550 / 3.825 Kč
MENICARE Pure 250 ml	9,00 EUR	229 Kč